

Preprocessing of fMRI data (basic)

Practical session
SPM Course 2016, Zurich

Andreea Diaconescu, Maya Schneebeli, Jakob Heinzle,
Lars Kasper, and Jakob Siemerikus

Translational Neuromodeling Unit (TNU)
Institute for Biomedical Engineering (IBT)
University and ETH Zürich



Translational Neuromodeling Unit





University of
Zurich^{UZH}

ETH zürich


Goals of this session

- Go through a preprocessing pipeline in SPM.
- Learn how to check whether some basic steps worked.
- Some basic file operations in SPM.
- Save, load and modify batches
- How to make your own preprocessing script.
- Answers to "all" your questions.


Preprocessing tools on the SPM GUI and Batch Editor



The Dataset: Event-related fMRI



- **Goal:** Investigate Repetition Suppression
- **How:** Each face presented twice during the session, 26 famous and 26 non-famous faces
 - 2x2 factorial design
 - Factor Fam(ilarity): long-term memory,
 - Level: Famous or Unfamous
 - Factor Rep(etition)
 - Level: 1 or 2
 - Task: Button press to decide fame




Stimulus Onsets

[all_conditions.mat](#)

*R. Henson et al.,
Cereb Cortex 2002*


Slice-Timing Correction (Temporal Preproc)

fMRI time-series



fmri.nii



Slice-Timing Correction




Slice-timing

afmri.nii

- **Goal:** Correct for different acquisition time of each slice within an image volume
- **How:** All voxel time series are aligned to acquisition time of 1 slice via Sinc-interpolation of each voxel's time series




Spatial Preprocessing



Realignment

fMRI time-series



REALIGN


- **Goal:** Correct for subject motion between volumes by minimising mean-squared difference
- **How:** Rigid-body transformation
 - Note: Realignment improves if images are *reoriented* in advance (find the origin, change header, use check reg with the canonical image)

fmri.nii

fmri.mat

rp_fmri.txt

meanfmri.nii



$$\begin{pmatrix} m_{11} & m_{12} & m_{13} & m_{14} \\ m_{21} & m_{22} & m_{23} & m_{24} \\ m_{31} & m_{32} & m_{33} & m_{34} \\ \mathbf{0} & \mathbf{0} & \mathbf{0} & \mathbf{1} \end{pmatrix}$$

Realignment parameters

(Headers changed)




Mean functional

Motion corrected




Co-registration



- Goal:** Match geometry of functional and structural images from same subject
 - Structural = high-resolution, geometrically correct image for later normalisation
- How:** Find affine transformation (rotation/translation/shear/scaling) that maximizes mutual information (similarity) between both images
 - Note: The role of functional and structural image could be reversed in this operation
 - Reference* = header unchanged
 - Source* = header will be changed


Normalisation: Unified Segmentation

Anatomical
MRI




struct.nii

Standard Space
Tissue Probability Maps



spm12/tpm/TPM.nii

Co-registered
functional




fmri.nii

**NORMALIZE
WRITE**

SEGMENT

c1-3struct.nii




Segmented

mstruct.nii




Bias-corrected
Structural

y_struct.nii



Deformation
Fields

wfmri.nii



Normalized Functional
(MNI Space)

wmstruct.nii




Normalized
Structural

Normalisation I: Compute Segmentation/Deformation Fields

Anatomical


MRI



struct.nii

Standard Space


Tissue Probability Maps



spm12/tpm/TPM.nii


SEGMENT

c1-3struct.nii




Segmented

mstruct.nii



Bias-corrected
Structural


y_struct.nii



Deformation
Fields

- **Goal:** Match geometry of subject brain to standard space (for group studies)
 - Integrated with segmentation, since mapping individual tissue classes is more robust
- **How:** Find non-linear transformation (deformation field) that makes tissue class distribution in structural image most plausible
 - Assuming coil inhomogeneity (bias), and deformations of reasonable anatomy
 - Note: uses a Bayesian maximum a posteriori estimation, where standard space tissue probability maps (TPMs) are the priors

Normalisation II: Warp via Write




- **Goal:** Write out functional/structural image in standard space for multi-subject statistical analysis

- to report findings in common anatomical space (e.g. MNI)


- **How:** Applies estimated deformation fields from Unified Segmentation to all functional and structural data

- Deformation fields have a displacement vector for each image voxel that tells where it should be moved to in standard space

Smoothing






Normalized Functional
(MNI Space)



SMOOTH



swfmri.nii



Smoothed Images

(Ready for Stats!)

- **Goal:** Increase sensitivity by reducing thermal noise and inter-subject variability in functional images
 - Plus a couple of more sophisticated reasons, e.g. matched-filter theorem, smoothness of residual random field for FWE-correction
- **How:** Convolution ("Blurring") with a 3D Gaussian kernel
 - Each voxel effectively becomes a weighted average of its surroundings (interpolation)
 - Typical full width at half maximum: 2-3 voxel

SPM Image File Prefixes `xfmri.nii/mat`

Prefix	Meaning	Typically applied to
a	Slice-timing corrected	functional
m	Bias-field corrected (<i>modulated</i>)	structural
mean	<i>Mean</i> of time-series	functional
r	<i>Re</i> -sliced (change image matrix, reset hdr)	both
s	Smoothed	functional
u	<i>Un</i> warped (during realign/FieldMap Tool)	functional
y_	Deformation field	structural
w	Normalized (<i>w</i> arped with deformation field)	both

SPM File Operations in Matlab Command Window

Operation	Matlab/SPM command
Retrieve multiple file names	<pre>fNames = spm_select('FPList', pwd, '^sM.*\.img\$');</pre>
Plot Image(s) with SPM	<pre>spm_check_registration(fNames)</pre>
Load image header(s)	<pre>hdr = spm_vol(fNames);</pre>
Load Image Matrix	<pre>X = spm_read_vols(hdr);</pre>
Inspect Batch in Batch Editor	<pre>spm_jobman('interactive', ... matlabbatch);</pre>
Run Batch from command line	<pre>spm_jobman('initcfg'); %1st time spm_jobman('run',matlabbatch);</pre>